





FEATURES


- Stylish diamond-dot pattern for everyday wear
- Durable construction with reinforced heel
- Excellent breathability and moisture management for added comfort

BEST SUITED FOR

- Patients who prefer fashionable options
- Pre-vein treatment
- Post-vein treatment
- Ideal for business or dress



20-30mmHg Closed Toe Pantyhose
Small Long in Graphite
712PSLW91



15-20mmHg Closed Toe Thigh-high
Medium Long in Black
711NMLW99

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Graphite (91)	Black (99)
THIGH W/GRIP-TOP (N)	●	●
PANTYHOSE (P)	●	●

COMPRESSION LEVELS

15-20mmHg	20-30mmHg
-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"-8.5" (18-21.5cm)	MEDIUM 8.5"-10" (21.5-25.5cm)	LARGE 10"-11.5" (25.5-29cm)			
CALF CIRCUMFERENCE	SMALL 11"-15.5" (28-39cm)	MEDIUM 12.5"-17.5" (31.5-45cm)	LARGE 14"-20" (35.5-51cm)			
THIGH CIRCUMFERENCE	SMALL 18-26" (45.5-66cm)	MEDIUM 19-28.5" (48-73cm)	LARGE 20.5-33" (52-84cm)			
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5-44.5" (83-112.5cm)	MEDIUM 34.5-50" (87-127cm)	LARGE 38.5-59.5" (98-150.5cm)			
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5-9.5, LONG LEG 8-12