



SEA ISLAND COTTON

for women & men


SEA ISLAND COTTON

FEATURES

- Knitted from Sea Island Cotton, the world's rarest cotton, for a fashionable look and a soft, luxurious feel
- Breathable and comfortable for sensitive skin
- Extended range of sizes to fit more patients

BEST SUITED FOR

- Patients who prefer natural fibers for everyday wear
- Management of venous insufficiency symptoms
- Long-distance travel



20-30mmHg Closed Toe Calf
Small Long in Brown
222CSLW11



20-30mmHg Closed Toe Calf
Medium Long in Navy
222CMLM10

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Navy(10)	Brown(11)	Black(99)
CALF (C)	●	●	●

COMPRESSION LEVELS

20-30mmHg
ALSO AVAILABLE IN 15-20mmHg, SEE PAGE 38

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"-8.5" (18-21.5cm)	MEDIUM 8.5"-10" (21.5-25.5cm)	LARGE 10"-11.5" (25.5-29cm)	EXTRA LARGE 11.5"-14" (29-36cm)				
CALF CIRCUMFERENCE	SMALL 11"-15.5" (28-39cm)	MEDIUM 12.5"-17.5" (31.5-45cm)	LARGE 14"-20" (35.5-51cm)	EXTRA LARGE 16"-24" (41-61cm)				
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7-11, LONG LEG 9.5-14 WOMEN'S SHOE SIZES: SHORT LEG 5.5-9.5, LONG LEG 8-12

ZURICH
COLLECTION

SIGVARIS MEDICAL

■ ALL-SEASON MERINO WOOL
for women & men




FEATURES

- Soft, comfortable, and itch-free
- Made from the highest quality extra-fine Australian Merino Wool
- Natural properties of wool make these socks odor resistant
- Thermoregulating for year-round wear and keeps feet dry

BEST SUITED FOR

- Patients who prefer fashionable and premium fabric options
- Management of venous insufficiency symptoms
- Long-distance travel



20–30mmHg Closed Toe Calf
Medium Long in Black
242CMLM99



20–30mmHg Closed Toe Calf
Small Long in Navy
242CSLW10

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Navy (10)	Brown (11)	Black (99)
CALF (C)	●	●	●

COMPRESSION LEVELS

20–30mmHg
ALSO AVAILABLE IN 15–20mmHg,
SEE PAGE 40

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)		EXTRA LARGE 11.5"–14" (29–36cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)		EXTRA LARGE 16"–24" (41–61cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 **WOMEN'S SHOE SIZES:** SHORT LEG 5.5–9.5, LONG LEG 8–12

ZURICH
COLLECTION

■ SIGVARIS MEDICAL