




COTTON
for women & men

FEATURES

- Only Supima Cotton touches the skin to ensure comfort
- Breathable and durable
- Classic ribbed look
- Ideal for everyday wear and for sensitive skin

BEST SUITED FOR

- Sensitive skin
- All-day cotton comfort in all types of climates
- Ideal for dress or casual shoes
- Ideal for long-distance travel



20–30mmHg Closed Toe Calf w/Grip-top
Medium Long in Crispa
232CMLW66/S



20–30mmHg Closed Toe Thigh-high
Medium Long in Black
232NMLM99

AVAILABLE STYLES AND COLORS

Closed Toe = ●	White (00)	Navy (10)	Black Mist (14)	Crispa (66)	Chocolate* (88)	Black (99)
Open Toe = ○						
CALF (C)	●	●	●	●○	●	●
CALF W/GRIP-TOP (C/S)				●*		●
THIGH W/GRIP-TOP (N)				●○		●

* 20–30mmHg ONLY

COMPRESSION LEVELS

20–30mmHg	30–40mmHg
-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)	MEDIUM 8.5"–10" (21.5–25.5cm)	LARGE 10"–11.5" (25.5–29cm)	EXTRA LARGE 11.5"–14" (29–36cm)				
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)	MEDIUM 12.5"–17.5" (31.5–45cm)	LARGE 14"–20" (35.5–51cm)	EXTRA LARGE 16"–24" (41–61cm)				
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18–26" (45.5–66cm)	MEDIUM 19–28.5" (48–73cm)	LARGE 20.5–33" (52–84cm)	XS & XL only available in Calf styles				
LEG LENGTH: WOMEN & OPEN TOE	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)					LONG >28.5" (72cm)
LEG LENGTH: MEN	SHORT ≤30" (76cm)	LONG >30" (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)	SHORT ≤30" (76cm)	LONG >30" (76cm)		
SIZE	SS	SL	MS	ML	LS	LL	XS	XL

MEN'S SHOE SIZES: SHORT LEG 7–11, LONG LEG 9.5–14 WOMEN'S SHOE SIZES: SHORT LEG 5.5–9.5, LONG LEG 8–12