




FEATURES

- Fashionable, sheer medical graduated compression hosiery
- Constructed from double-covered inlay yarns for softness against the skin
- Excellent breathability and moisture management for added comfort

BEST SUITED FOR

- Patients who prefer fashionable options
- Pre-vein treatment
- Post-vein treatment
- Ideal for business or dress



15–20mmHg Closed Toe Pantyhose
Small Long in Black
781PSLW99



20–30mmHg Closed Toe Thigh-high
Medium Long in Suntan
782NMLW36

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Dark Navy (08)	Natural (33)	Suntan (36)	Mocha (85)	Nightshade (94)	Black (99)
Open Toe = ○						
CALF (C)	●	●●	●●	●	●	●●
THIGH W/GRIP-TOP (N)	●	●●	●●	●	●	●●
PANTYHOSE (P)	●	●●	●●	●	●	●●

COMPRESSION LEVELS

15–20mmHg	20–30mmHg	30–40mmHg
-----------	-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"–8.5" (18–21.5cm)		MEDIUM 8.5"–10" (21.5–25.5cm)		LARGE 10"–11.5" (25.5–29cm)	
CALF CIRCUMFERENCE	SMALL 11"–15.5" (28–39cm)		MEDIUM 12.5"–17.5" (31.5–45cm)		LARGE 14"–20" (35.5–51cm)	
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18–26" (45.5–66cm)		MEDIUM 19–28.5" (48–73cm)		LARGE 20.5–33" (52–84cm)	
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5–44.5" (83–112.5cm)		MEDIUM 34.5–50" (87–127cm)		LARGE 38.5–59.5" (98–150.5cm)	
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5–9.5, LONG LEG 8–12