

NATURAL RUBBER
for women & men

NATURAL RUBBER



FEATURES

- Fits a broad range of medical patients, including those with extra large ankle and calf circumferences
- Therapeutic solution for those with long-term therapy needs
- Minimal fiber fatigue throughout the day
- This product is often combined with our Cotton Underliner (see Accessories pages) for patient comfort

BEST SUITED FOR

- Severe varicosities
- Severe edema
- Stasis
- Post-ulcer treatment
- Lymphedema
- Venous ulcers
- C3-C6 (CEAP)



30-40mmHg Open Toe Thigh-high w/Grip-top
Small Average Long in Beige
503NS2077



40-50mmHg Open Toe Calf
Medium Full Long in Beige
504CM4077

AVAILABLE STYLES AND COLOR

Open Toe = <input checked="" type="radio"/>	BEIGE (77)	
CALF (C)	<input type="radio"/>	CALF ONLY AVAILABLE IN S, M, AND L FULL CIRCUMFERENCE FOR 50-60mmHg
THIGH (T) WITHOUT GRIP-TOP	<input type="radio"/>	THIGH ALSO AVAILABLE IN XS
THIGH (N) WITH GRIP-TOP	<input type="radio"/>	
PANTYHOSE (P)	<input type="radio"/>	PANTYHOSE ONLY AVAILABLE IN 30-40mmHg AND ONLY IN FULL THIGH CIRCUMFERENCE
THIGH (W) WITH WAIST ATTACHMENT (INDICATE LEG WITH /L OR/R)	<input type="radio"/>	THIGH WITH WAIST ATTACHMENT ONLY AVAILABLE IN 30-40mmHg AND 40-50mmHg ONLY AVAILABLE IN FULL THIGH CIRCUMFERENCE

*CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS. FOR ALL INDICATIONS AND CONTRAINDICATIONS, PLEASE SEE SIGVARIS PRESCRIPTION PAD. SIGVARIS ACCEPTS NO LIABILITY FOR NON-OBSERVANCE OF CONTRAINDICATIONS OR CAUTIONS.

COMPRESSION LEVELS

30-40mmHg	40-50mmHg	50-60mmHg
-----------	-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL				MEDIUM				LARGE				X-LARGE (CALF ONLY)			
	7.5"-9" (19-23cm)				9"-12" (23-30cm)				12"-15.5" (30-40cm)				15.5"-18" (40-46cm)			
CALF CIRCUMFERENCE	AVERAGE		FULL		AVERAGE		FULL		AVERAGE		FULL		AVERAGE		FULL	
	11"-15" (28-38cm)	≤17" (≤43cm)	12.5"-17" (32-43cm)	≤19" (≤48cm)	14.5"-19" (36-48cm)	≤22.5" (≤57cm)	18"-22.5" (45-57cm)	≤27" (≤68cm)								
CALF CIRCUMFERENCE FOR 505—FULL ONLY	505C: FULL CALF ONLY		12.5"-17" (32-43cm)		505C: FULL CALF ONLY		14.5"-19" (36-48cm)		505C: FULL CALF ONLY		18"-22.5" (45-57cm)		505C NOT AVAILABLE IN XL		505C NOT AVAILABLE IN XL	
	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG
CALF LENGTH	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)	≤15.5" (40cm)	>15.5" (40cm)
	AVERAGE		FULL		AVERAGE		FULL		AVERAGE		FULL		CUSTOM SIZES AVAILABLE. PLEASE CONTACT CUSTOMER CARE FOR DETAILS			
THIGH CIRCUMFERENCE	16"-21" (41-54cm)		≤24" (≤61cm)		19"-24" (48-61cm)		≤28" (≤71cm)		23"-28" (58-71cm)		≤32" (≤81cm)					
	505T: FULL THIGH ONLY		≤24" (≤61cm)		505T: FULL THIGH ONLY		≤28" (≤71cm)		505T: FULL THIGH ONLY		≤32" (≤81cm)					
LEG LENGTH	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	SHORT	LONG	CUSTOM SIZES AVAILABLE. PLEASE CONTACT CUSTOMER CARE FOR DETAILS			
	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)	≤28" (71cm)	>28" (71cm)				
SIZE	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4				

GUIDELINES FOR HIP CIRCUMFERENCES: SMALL AVERAGE <42" (<106.5cm), SMALL FULL <48" (<122cm), MEDIUM AVERAGE <48" (<122cm), MEDIUM FULL <54" (<137cm), LARGE AVERAGE <55" (<140cm), LARGE FULL <61" (<155cm)