



SELECT COMFORT


for women closed toe and men & women open toe

FEATURES

- Wide range of sizing options to fit every body type
- Petite and Plus sizing available
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material allows for ease in putting on and taking off

BEST SUITED FOR

- Post-vein treatment (stripping, ligation, phlebectomy, sclero, EVLT, RF ablation)
- Pregnancy
- Options for petite and plus-sized patients



SELECT COMFORT



30-40mmHg Closed Toe Thigh-high Medium Long in Black
863NMLW99



20-30mmHg Open Toe Calf with Grip-top Small Short in Crispa
862CSS066/S

AVAILABLE STYLES AND COLORS

	Natural (33)	Suntan (36)	Crispa (66)	Black (99)
Closed Toe = ●				
Open Toe = ○				
CALF (C)	●	●	○	●●
CALF W/GRIP-TOP (C/S)			●●	●
THIGH W/GRIP-TOP (N)	●	●	○	●
PANTYHOSE (P)	●	●		●
PANTYHOSE PLUS (P/P)			●	●
THIGH WITH WAIST ATTACHMENT (W) (INDICATE LEG WITH /L OR /R)			○	

NOTE: 860 OPEN TOE PRODUCTS ARE UNISEX

COMPRESSION LEVELS

20-30mmHg	30-40mmHg
-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"-8.5" (18-21.5cm)			MEDIUM 8.5"-10" (21.5-25.5cm)			LARGE 10"-11.5" (25.5-29cm)		EXTRA LARGE 11.5"-14" (29-36cm)	
CALF CIRCUMFERENCE	SMALL 11"-15.5" (28-39cm)			MEDIUM 12.5"-17.5" (31.5-45cm)			LARGE 14"-20" (35.5-51cm)		EXTRA LARGE 16"-24" (41-61cm)	
CALF LENGTH	PETITE <14" (35.5cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	PETITE <14" (35.5cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18-26" (45.5-66cm)			MEDIUM 19-28.5" (48-73cm)			LARGE 20.5-33" (52-84cm)		XS & XL only available in Calf styles	
LEG LENGTH: WOMEN & OPEN TOE	PETITE ≤25.5" (65cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	PETITE ≤25.5" (65cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)		
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5-44.5" (83-112.5cm)			MEDIUM 34.5-50" (87-127cm)			LARGE 38.5-59.5" (98-150.5cm)			
SIZE	SP*	SS	SL	MP*	MS	ML	LS	LL	XS	XL

WOMEN'S SHOE SIZES: PETITE LEG 4-8, SHORT LEG 5.5-9.5, LONG LEG 8-12

PANTYHOSE PLUS SIZES HIP CIRCUMFERENCE GUIDELINES: SMALL: <56" (142cm), MEDIUM: <62" (157.5cm), LARGE: <70" (178cm)

*SP & MP SIZES ARE AVAILABLE IN CLOSED TOE ONLY, WITHOUT GRIP-TOP, AND ARE UNAVAILABLE IN PANTYHOSE PLUS SIZES