

SOFT OPAQUE
for women




FEATURES


- Comfortable tights with lasting microfiber softness
- Advanced moisture management keeps legs cool and dry
- Comfortable wide top band prevents knee-highs from rolling down

BEST SUITED FOR

- Patients who prefer fashionable options
- Post-vein treatment
- Pregnancy
- Long-distance travel



20-30mmHg Closed Toe Calf
Medium Long in Graphite
842CMLW91



15-20mmHg Closed Toe Thigh-high
Medium Long in Black
841NMLW99

AVAILABLE STYLES AND COLORS

	Midnight Blue (09)	Nude (35)	Pecan** (40)	Mulberry** (60)	Espresso (89)	Graphite (91)	Black (99)
Closed Toe = ●							
Open Toe = ○							
CALF (C)	●	●●	●	●	●	●	●●
THIGH W/GRIP-TOP (N)	●	●●	●	●	●	●	●●
PANTYHOSE (P)	●	●●	●	●	●	●	●●
MATERNITY PANTYHOSE (M)*							●

* 20-30mmHg ONLY
** 15-20 & 20-30mmHg ONLY

COMPRESSION LEVELS

15-20mmHg	20-30mmHg	30-40mmHg
-----------	-----------	-----------

SIZING CHART

ANKLE CIRCUMFERENCE	SMALL 7"-8.5" (18-21.5cm)	MEDIUM 8.5"-10" (21.5-25.5cm)	LARGE 10"-11.5" (25.5-29cm)			
CALF CIRCUMFERENCE	SMALL 11"-15.5" (28-39cm)	MEDIUM 12.5"-17.5" (31.5-45cm)	LARGE 14"-20" (35.5-51cm)			
CALF LENGTH	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)	SHORT ≤16" (41cm)	LONG >16" (41cm)
THIGH CIRCUMFERENCE	SMALL 18-26" (45.5-66cm)	MEDIUM 19-28.5" (48-73cm)	LARGE 20.5-33" (52-84cm)			
LEG LENGTH: WOMEN	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)	SHORT ≤28.5" (72cm)	LONG >28.5" (72cm)
GUIDELINES FOR HIP CIRCUMFERENCE	SMALL 32.5-44.5" (83-112.5cm)	MEDIUM 34.5-50" (87-127cm)	LARGE 38.5-59.5" (98-150.5cm)			
SIZE	SS	SL	MS	ML	LS	LL

WOMEN'S SHOE SIZES: SHORT LEG 5.5-9.5, LONG LEG 8-12

MATERNITY HIP CIRCUMFERENCE GUIDELINES: SMALL: <56" (142cm), MEDIUM: <62" (157.5cm), LARGE: <70" (178cm)