




CASUAL COTTON
for women & men

FEATURES

- Made with premium Supima® Cotton for optimum softness and durability
- Flat, low friction toe seam to prevent irritation
- Soft, stretchable knit-in band prevents pinching

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Preventative (travel, occupational risk, pregnancy, family history and comfort)



146: WOMEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size B in Khaki
146CB30

186: MEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size C in Black
186CC99

AVAILABLE STYLE AND COLORS

Closed Toe = ●	White (00)	Navy (10)	Brown (11)	Khaki (30)	Black (99)
CALF (C)	●	●	●	●	●
COTTON MATERNITY SOCKS (C/M)	●				●

COMPRESSION LEVEL

15–20mmHg

SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" (18–23cm)	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
SIZE	A	B	C