

CUSHIONED COTTON
for women & men




FEATURES

- Thick, padded sole, heel and toe to protect feet
- Flat, low friction toe seam
- Soft, stretchable knit-in band prevents pinching
- Ideal for everyday wear and for sensitive skin
- Suitable for patients with diabetes

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Preventative (travel, occupational risk, pregnancy, family history and comfort)



142: WOMEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size B in White
142CB00

182: MEN'S CALF STYLE



15–20mmHg Closed Toe Calf
Size C in Black
182CC99

AVAILABLE STYLE AND COLORS

Closed Toe = ●	White (00)	Black (99)
	CALF (C) ●	●

COMPRESSION LEVELS

15–20mmHg

ALSO AVAILABLE IN 20–30mmHg,
SEE PAGE 8

SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7"–9" (18–23cm)	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)
MEN'S SHOE	6.5–8.5	9–11	11.5–14
ANKLE CIRCUMFERENCE	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
SIZE	A	B	C