




FEATURES

- New updated look features a sheer toe and heel area
- Improved for softer feel against the skin
- Incredible sheer design in a variety of fashionable colors
- Comfortable stretch is excellent for everyday wear

BEST SUITED FOR

- Tired, heavy, restless, aching legs
- Preventative (travel, occupational risk, pregnancy, family history and comfort)



15–20mmHg Closed Toe Calf
Size B in Natural
120CB33



15–20mmHg Closed Toe Pantyhose
Size C in Black
120PC99

AVAILABLE STYLES AND COLORS

Closed Toe = ●	Charcoal (12)	Taupe (29)	Natural (33)	Suntan (36)	Café (73)	Black (99)
Open Toe = ○						
CALF (C)	●	●	●	●	●	●
THIGH (N)	●	●	●	●	●	●
PANTYHOSE (P)	●	●	●	●	●	●
MATERNITY PANTYHOSE (M)			●			●

COMPRESSION LEVEL

15–20mmHg



Silicone band grip-top on Thigh-high

CALF AND THIGH SIZING CHART

WOMEN'S SHOE	5–7	7.5–9.5	10–12
ANKLE CIRCUMFERENCE	7.5"–9" (19–23cm)	8"–10" (20.5–25cm)	8.5"–10.5" (22–27cm)
SIZE	A	B	C

WEIGHT GUIDELINES FOR THIGH STYLE:
SIZE A: <135LBS., SIZE B: <180LBS., SIZE C: <230LBS.

PANTYHOSE AND MATERNITY SIZING CHART

	100/45	105/48	110/50	115/52	120/54	125/57	130/59	135/61	140/64	145/66	150/68	155/70	160/73	165/75	170/77	175/80	180/82	185/84	190/86	195/89	200/91	205/93	210/95	215/98	220/100	225/102	230/105	
4'9" / 145cm																												
4'10" / 147cm																												
4'11" / 150cm																												
5'0" / 152cm																												
5'1" / 155cm																												
5'2" / 157cm																												
5'3" / 160cm																												
5'4" / 162cm																												
5'5" / 165cm																												
5'6" / 167cm																												
5'7" / 170cm																												
5'8" / 173cm																												
5'9" / 175cm																												
5'10" / 178cm																												
5'11" / 180cm																												
6' / 183cm																												

NOTE: PLEASE USE PRE-PREGNANCY WEIGHT FOR MATERNITY PANTYHOSE.