

Juzo® Basic Series



Economically Priced



Two-way stretch elasticity



Available knee-high, thigh-high and pantyhose



Freedom in Motion™


Juzo® Basic Series



Juzo® Basic


An economically priced compression garment

Model	4410, 4411, 4412
Compression	15-20, 20-30, 30-40 mmHg
Styles	knee-high (AD), thigh-high (AG), pantyhose (AT)
Regular Sizes	I - IV, regular and short length
Foot	open toe or full foot
AG Silicone Border	beaded

Colors Available
 beige

Size Chart For Juzo® Basic Compression Stockings and Pantyhose

4410, 4411 & 4412 15-20, 20-30 and 30-40 mmHg



	I	II	III	IV
g	44-66 cm 17 1/4"-26"	48-71 cm 18 3/4"-28"	52-76 cm 20 1/2"-29 3/4"	56-81 cm 22"-32"
g with silicone	44-62 cm 17 1/4"-24 1/2"	48-67 cm 18 3/4"-26 1/4"	52-72 cm 20 1/2"-28 1/4"	56-77 cm 22"-30 1/4"
c	28-38 cm 11"-15"	31-42 cm 12"-16 1/2"	34-46 cm 13 1/4"-18"	37-50 cm 14 1/2"-19 3/4"
b	18-21.5 cm 7"-8 1/2"	21.5-25 cm 8 1/2"-9 3/4"	25-28 cm 9 3/4"-11"	28-31.5 cm 11"-12 1/2"

length A-d
A-D Regular - 15 3/4"-18 1/4" / **40-46 cm**
A-D Short - 13-15 3/4" / **34-40 cm**

length A-g
A-G Regular and A-T Regular - 28 1/4"-32 3/4" / **72-83 cm**
A-G Short and A-T Short - 24 3/4"-28 1/4" / **63-72 cm**

For more information
www.juzousa.com
 or call 1-800-222-4999

