

# JOBST® Opaque

**NEW!**  
Midnight Navy



**Softness and ease  
of donning for  
exceptionally  
comfortable leg  
therapy.**

**Medical Compression Stockings  
15-20mmHg\*, 20-30mmHg\*, 30-40mmHg\***

#1 physician recommended

**JOBST®**

Comfort, Health and Style!



## Features

## Benefits



Enhanced with a silicone blend softener

Unmatched softness and superior moisture wicking for cool and comfortable legs.

Contoured panty and anatomical foot design with a relaxed toe cap and enlarged heel pocket

Opaque is easier to don than any other stocking you have worn before.

Refined opacity

Ensures a perfect fit from hip to toe all day long for exceptional wearing comfort.

Selected, strong yarns and a reinforced heel and toe


Covers blemishes and scars giving your leg an elegant appearance.

Fashionable but durable backed by Our Total Satisfaction Guarantee.

## What is JOBST® Advanced Comfort

JOBST® Advanced Comfort combines exceptional yarns with a unique finish to create unbelievable wearing comfort. This technology encompasses:

- Unmatched Softness: finest yarns and unique finishes are combined to give a perfect soft touch
- Superior Moisture Wicking: Jobst Advanced Comfort boosts the spreading speed of any moisture on the fabric so that it evaporates faster. This keeps your legs comfortable all day long.



#1 physician recommended



Comfort, Health and Style!



*Want softer compression therapy? Here you go ...*


## **Why JOBST® Opaque?**

- JOBST® Opaque is the softest Opaque ever and JOBST® Advanced Comfort Opaque gives unbelievable wearing comfort.
- The contoured panty provides a perfect fit around the hips to make our waist highs even more comfortable.
- JOBST®'s silicone blend softener makes Opaque very easy to don.
- New elegant collection of colors to mix and match with all wardrobes of any season.

## **Women voted. We delivered.**

Women have voted! Trials have proven that women clearly prefer the new Opaque. 90% of women prefer the new Opaque vs. the current Opaque due to its perfect soft touch\*.

\* Blinded preference trial carried out by JOBST®.  
Data available upon request.



**A variety of colors and styles to choose from.**

## Colors



*Classic Black*

*Midnight Navy*

*Espresso*

*Honey*

*Natural*

Actual garment color may vary from image shown.

## Styles



*Knee High  
Closed Toe*



*Knee High  
Open Toe*



*Thigh High  
Closed Toe*



*Thigh High  
Open Toe*



*Pantyhose*

## Knee High ordering information:

**Knee High  
(Closed Toe)**



Size	Catalog Number/Color					Petite (Closed Toe)	Size	Catalog Number/Color		Compression mmHg*
	Classic Black	Honey	Natural	Espresso	<b>NEW!</b> Midnight Navy			Classic Black	Natural	
S	115200	115681	115212	115685	115752		S	115600	115604	15-20 mmHg
M	115201	115682	115213	115686	115753		M	115601	115605	15-20 mmHg
L	115202	115683	115214	115687	115754		L	115602	115606	15-20 mmHg
XL	115203	115684	115215	115688	115755		XL	115603	115607	15-20 mmHg
S	115132	115456	115270	115460	115740		S	115608	115614	20-30 mmHg
M	115133	115457	115271	115461	115741		M	115609	115615	20-30 mmHg
L	115134	115458	115272	115462	115742		L	115610	115616	20-30 mmHg
XL	115135	115459	115273	115463	115743		XL	115611	115617	20-30 mmHg
LFC	115364		115378							20-30 mmHg
XLFC	115365		115379							20-30 mmHg
S	115168	115701	115282	115705			S	115626	115632	30-40 mmHg
M	115169	115702	115283	115706			M	115627	115633	30-40 mmHg
L	115170	115703	115284	115707			L	115628	115634	30-40 mmHg
XL	115171	115704	115285	115708			XL	115629	115635	30-40 mmHg
LFC	115368		115370							30-40 mmHg
XLFC	115369		115371							30-40 mmHg

#1 physician recommended



Comfort, Health and Style!

This panel size is 7 3/4" wide to allow for hole punch at the fold in



### Knee High ordering information:

Knee High  
(Open Toe)



Catalog Number/Color			Catalog Number/Color Compression mmHg*		
Size	Classic Black	Natural	Petite (Open Toe) Size	Natural	Compression mmHg*
S	115335	115331	S	115069	15-20 mmHg
M	115336	115332	M	115070	15-20 mmHg
L	115337	115333	L	115071	15-20 mmHg
XL	115338	115334	XL	115072	15-20 mmHg
S	115384	115480	S	115620	20-30 mmHg
M	115385	115481	M	115621	20-30 mmHg
L	115386	115482	L	115622	20-30 mmHg
XL	115387	115483	XL	115623	20-30 mmHg
LFC	115376	115366			20-30 mmHg
XLFC	115377	115367			20-30 mmHg
S	115388	115496	S	115638	30-40 mmHg
M	115389	115497	M	115639	30-40 mmHg
L	115390	115498	L	115640	30-40 mmHg
XL	115391	115499	XL	115641	30-40 mmHg
LFC	115380	115382			30-40 mmHg
XLFC	115381	115383			30-40 mmHg

### Fiber contents:

15-20 mmHg\*: 78% Nylon and 22% Spandex/Elasthan; 20-30 mmHg\*: 78% Nylon and 22% Spandex/Elasthan; 30-40 mmHg\*: 70% Nylon and 30% Spandex/Elasthan

### Thigh High ordering information:

Thigh High  
(Closed Toe)  
(Silicone Dot Band)



Catalog Number/Color					<b>NEW!</b>	Catalog Number/Color Compression mmHg*		
Size	Classic Black	Honey	Natural	Espresso	Midnight Navy	Petite (Closed Toe) (Silicone Dot Band) Size	Natural	Compression mmHg*
S	115504	115689	115508	115693	115760	S	115652	15-20 mmHg
M	115505	115690	115509	115694	115761	M	115653	15-20 mmHg
L	115506	115691	115510	115695	115762	L	115654	15-20 mmHg
XL	115507	115692	115511	115696	115763	XL	115655	15-20 mmHg
S	115144	115464	115274	115468	115748	S	115660	20-30 mmHg
M	115145	115465	115275	115469	115749	M	115661	20-30 mmHg
L	115146	115466	115276	115470	115750	L	115662	20-30 mmHg
XL	115147	115467	115277	115471	115751	XL	115663	20-30 mmHg
S	115180	115709	115286	115713		S	115672	30-40 mmHg
M	115181	115710	115287	115714		M	115673	30-40 mmHg
L	115182	115711	115288	115715		L	115674	30-40 mmHg
XL	115183	115712	115289	115716		XL	115675	30-40 mmHg

### Thigh High ordering information:

Thigh High  
(Open Toe)  
(Silicone Dot Band)



Catalog Number/Color			Catalog Number/Color Compression mmHg*			
Size	Classic Black	Honey	Natural	Petite (Open Toe) (Silicone Dot Band) Size	Natural	Compression mmHg*
S	115556	115697	115552	S	115644	15-20 mmHg
M	115557	115698	115553	M	115645	15-20 mmHg
L	115558	115699	115554	L	115646	15-20 mmHg
XL	115559	115700	115555	XL	115647	15-20 mmHg
S	115560	115472	115544	S	115664	20-30 mmHg
M	115561	115473	115545	M	115665	20-30 mmHg
L	115562	115474	115546	L	115666	20-30 mmHg
XL	115563	115475	115547	XL	115667	20-30 mmHg
S	115564	115717	115548	S	115676	30-40 mmHg
M	115565	115718	115549	M	115677	30-40 mmHg
L	115566	115719	115550	L	115678	30-40 mmHg
XL	115567	115720	115551	XL	115679	30-40 mmHg


### Fiber contents:

15-20 mmHg\*: 72% Nylon, 18% Spandex, 10% Silicone; 20-30 mmHg\*: 70% Nylon, 16% Spandex, 14% Silicone; 30-40 mmHg\*: 62% Nylon, 28% Spandex, 10% Silicone

# JOBST® Opaque

15-20 mmHg\*, 20-30 mmHg\*, 30-40 mmHg\*







## Pantyhose ordering information:

Pantyhose	Catalog Number/Color			Compression mmHg*	
	Size	Classic Black	Natural		NEW! Midnight Navy
	S	115220	115232	115756	15-20 mmHg
	M	115221	115233	115757	15-20 mmHg
	L	115222	115234	115758	15-20 mmHg
	XL	115223	115235	115759	15-20 mmHg
	S	115156	115278	115744	20-30 mmHg
	M	115157	115279	115745	20-30 mmHg
	L	115158	115280	115746	20-30 mmHg
	XL	115159	115281	115747	20-30 mmHg
	S	115192	115290		30-40mmHg
	M	115193	115291		30-40mmHg
	L	115194	115292		30-40mmHg
	XL	115195	115293		30-40mmHg

## Fiber contents:

15-20 mmHg\*: 76% Nylon, 24% Spandex; 20-30 mmHg\*: 76% Nylon, 24% Spandex; 30-40 mmHg\*: 66% Nylon, 34% Spandex

## Sizing chart:

SIZE	ANKLE 	CALF 	THIGH 	HIP 
<b>S</b>	7"-8 <sup>1</sup> / <sub>4</sub> " (18-21 cm)	11"-15" (28-38 cm)	15 <sup>3</sup> / <sub>4</sub> "-24 <sup>3</sup> / <sub>8</sub> " (40-62 cm)	28"-46" (71-117 cm)
<b>M</b>	8 <sup>3</sup> / <sub>8</sub> "-9 <sup>7</sup> / <sub>8</sub> " (21-25 cm)	11 <sup>7</sup> / <sub>8</sub> "-16 <sup>1</sup> / <sub>2</sub> " (30-42 cm)	18 <sup>1</sup> / <sub>8</sub> "-27 <sup>1</sup> / <sub>2</sub> " (46-70 cm)	30"-50" (76-127 cm)
<b>L</b>	10"-11 <sup>3</sup> / <sub>8</sub> " (25-29 cm)	12 <sup>1</sup> / <sub>2</sub> "-18 <sup>1</sup> / <sub>8</sub> " (32-46 cm)	21 <sup>1</sup> / <sub>4</sub> "-30 <sup>3</sup> / <sub>4</sub> " (54-78 cm)	32"-54" (81-137 cm)
<b>XL</b>	11 <sup>1</sup> / <sub>2</sub> "-13" (29-33 cm)	13 <sup>3</sup> / <sub>8</sub> "-19 <sup>5</sup> / <sub>8</sub> " (34-50 cm)	23 <sup>5</sup> / <sub>8</sub> "-32" (60-81 cm)	40"-65" (102-166 cm)
<b>LFC</b>	10"-11 <sup>7</sup> / <sub>8</sub> " (25-30 cm)	18"-24" (46-61 cm)		
<b>XLFC</b>	12"-14" (30-36 cm)	18"-24" (46-61 cm)		
	Knee High Regular: Greater than 15" (38 cm) Knee High Petite: 15" (38 cm) or less			Thigh High Regular: Greater than 27" (68.5cm) Thigh High Petite: 27" (68.5cm) or less

For Knee High, Measure Ankle & Calf. For Thigh High, Measure Ankle, Calf & Thigh. For Waist High, Measure Ankle, Thigh & Hip.

## Care instructions:

Machine washable in lukewarm water. Do not bleach, tumble dry, dry clean, or iron.



15-20 mmHg\*  
Compression



20-30 mmHg\*  
Compression



30-40 mmHg\*  
Compression

\* The mean compression of an average ankle size.  
Please refer to the product label and /or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.jobst.com  
To order toll-free: JOBST (+1) 800 537 1063

60419 R2

© 2012 BSN medical Inc.

physician recommended



Comfort, Health and Style!

REV 03/12