

Relief®

15-20 mmHg*, 20-30 mmHg*, 30-40 mmHg*

Waist High ordering information:

Waist High
(Closed Toe)



Size	Catalog Number		Compression
	Beige	Beige	
S	114656	114656	20-30 mmHg
M	114657	114657	20-30 mmHg
L	114658	114658	20-30 mmHg
XL	114659	114659	20-30 mmHg
S	114664	114664	30-40 mmHg
M	114665	114665	30-40 mmHg
L	114666	114666	30-40 mmHg
XL	114667	114667	30-40 mmHg

Waist High
(Open Toe)



Size	Catalog Number		Compression
	Beige	Beige	
S	114660	114660	20-30 mmHg
M	114661	114661	20-30 mmHg
L	114662	114662	20-30 mmHg
XL	114663	114663	20-30 mmHg
S	114668	114668	30-40 mmHg
M	114669	114669	30-40 mmHg
L	114670	114670	30-40 mmHg
XL	114671	114671	30-40 mmHg

Fiber contents:

20-30 mmHg & 30-40 mmHg: 75% Nylon, 25% Spandex, WH CT
20-30 mmHg & 30-40 mmHg: 78% Nylon, 22% Spandex, WH OT

Chaps ordering information:

Chaps - Left Leg
(Open Toe)



Size	Catalog Number		Compression
	Beige	Beige	
S	114672	114672	20-30 mmHg
M	114673	114673	20-30 mmHg
L	114674	114674	20-30 mmHg
XL	114675	114675	20-30 mmHg
S	114784	114784	30-40 mmHg
M	114785	114785	30-40 mmHg
L	114786	114786	30-40 mmHg
XL	114787	114787	30-40 mmHg

Chaps - Right Leg
(Open Toe)



Size	Catalog Number		Compression
	Beige	Beige	
S	114676	114676	20-30 mmHg
M	114677	114677	20-30 mmHg
L	114678	114678	20-30 mmHg
XL	114679	114679	20-30 mmHg
S	114788	114788	30-40 mmHg
M	114789	114789	30-40 mmHg
L	114790	114790	30-40 mmHg
XL	114791	114791	30-40 mmHg

Fiber contents:

20-30 mmHg & 30-40 mmHg: 78% Nylon, 22% Spandex, Chaps OT

Chaps - Both Legs
(Open Toe)



Size	Catalog Number		Compression
	Beige	Beige	
S	114680	114680	20-30 mmHg
M	114681	114681	20-30 mmHg
L	114682	114682	20-30 mmHg
XL	114683	114683	20-30 mmHg
S	114792	114792	30-40 mmHg
M	114793	114793	30-40 mmHg
L	114794	114794	30-40 mmHg
XL	114795	114795	30-40 mmHg



* The mean compression for an average size ankle.

NOTE: Please refer to the product label and/or package insert for full instructions on the safe use of these products.



BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.jobst.com
To order toll-free: JOBST (+1) 800 537 1063

52211 R8

© 2009 BSN medical Inc.

Printed in U.S.A.



Comfort, Health and Style!

REV 05/09

New!
Knee High
with Silicone Band

Relief®

- Leg Therapy
- Comfortable
- Durable
- Compliance

Medical Compression Stockings
15-20 mmHg*, 20-30 mmHg*, 30-40 mmHg*



Comfort, Health and Style!

Relief®

15-20 mmHg*, 20-30 mmHg*, 30-40 mmHg*

Quality and therapy at a moderate price

JOBST® Relief delivers compression therapy at an affordable price. This is especially important when reimbursement, provided by managed care, limits product choices.

Features & benefits:

Leg Therapy

- Accurate gradient compression therapy provides relief from tired, aching legs and varicose veins.

Comfortable

- 3-D knit structure uses air-covered spandex yarn that is soft to touch and comfortable for all-day wear.
- Closed toe styles are designed with a roomy toe for patients who require additional “wiggle room”.

Durable

- Reinforced heel is designed for long lasting wear.

Compliance

- Knee, thigh and waist high styles available in closed or open toe.
- Relief Knee High also comes in a full calf size.
- Chap style available in open toe.

Available colors:

Beige








Black



Actual garment color may vary from printed swatch sample.

Sizing chart:

SIZE	ANKLE 	CALF 	THIGH 	HIP 
S	7"-8 ¹ / ₄ " (18-21 cm)	11"-15" (28-38 cm)	15 ³ / ₄ "-24 ³ / ₈ " (40-62 cm)	28"-46" (71-117 cm)
M	8 ³ / ₈ "-9 ⁷ / ₈ " (21-25 cm)	11 ⁷ / ₈ "-16 ¹ / ₂ " (30-42 cm)	18 ¹ / ₈ "-27 ¹ / ₂ " (46-70 cm)	30"-50" (76-127 cm)
L	10"-11 ³ / ₈ " (25-29 cm)	12 ¹ / ₂ "-18 ¹ / ₈ " (32-46 cm)	21 ¹ / ₄ "-30 ³ / ₄ " (54-78 cm)	32"-54" (81-137 cm)
XL	11 ¹ / ₂ "-13" (29-33 cm)	13 ³ / ₈ "-19 ⁵ / ₈ " (34-50 cm)	23 ⁵ / ₈ "-32" (60-81 cm)	40"-65" (102-165 cm)
LFC	10"-11 ⁷ / ₈ " (25-30 cm)	18"-24" (46-61 cm)		
XLFC	12"-14" (30-36 cm)	18"-24" (46-61 cm)		
	Regular: Greater than 15" (38cm). Available in all above sizes.			
	Petite: 15" (38cm) or less. *M & L available in petite sizes.			

Knee High: Measure ankle and calf. Thigh High: Measure ankle, calf and thigh. Waist High and Chaps: Measure ankle, thigh and hip.
LFC: Large Full Calf XLFC: Extra-Large Full Calf.

Knee High ordering information:

Knee High with Silicone (Closed Toe)



Size	Catalog Number	Compression
M	114743	20-30 mmHg
L	114744	20-30 mmHg
XL	114745	20-30 mmHg
LFC	114746	20-30 mmHg
XLFC	114747	20-30 mmHg

Knee High with Silicone (Open Toe)



Size	Catalog Number	Compression
M	114748	20-30 mmHg
M-Petite	114760	20-30 mmHg
L	114749	20-30 mmHg
L-Petite	114761	20-30 mmHg
XL	114750	20-30 mmHg
LFC	114751	20-30 mmHg
XLFC	114752	20-30 mmHg
M	114753	30-40 mmHg
L	114754	30-40 mmHg
XL	114755	30-40 mmHg
LFC	114756	30-40 mmHg
XLFC	114757	30-40 mmHg

Knee High ordering information:

Knee High (Closed Toe)



Size	Catalog Number		Compression
	Beige	Black	
S	114806	114812	15-20 mmHg
M	114807	114813	15-20 mmHg
L	114808	114814	15-20 mmHg
XL	114809	114815	15-20 mmHg
LFC	114810	114816	15-20 mmHg
XLFC	114811	114817	15-20 mmHg
S	114620	114730	20-30 mmHg
M	114621	114731	20-30 mmHg
L	114622	114732	20-30 mmHg
XL	114623	114733	20-30 mmHg
LFC	114698	114734	20-30 mmHg
XLFC	114624	114735	20-30 mmHg
S	114630	114736	30-40 mmHg
M	114631	114737	30-40 mmHg
L	114632	114738	30-40 mmHg
XL	114633	114739	30-40 mmHg
LFC	114699	114740	30-40 mmHg
XLFC	114634	114741	30-40 mmHg

Knee High (Open Toe)



Size	Catalog Number		Compression
	Beige	Black	
S	114800		15-20 mmHg
M	114801		15-20 mmHg
L	114802		15-20 mmHg
XL	114803		15-20 mmHg
LFC	114804		15-20 mmHg
XLFC	114805		15-20 mmHg
S	114625		20-30 mmHg
M	114626		20-30 mmHg
L	114627		20-30 mmHg
XL	114628		20-30 mmHg
LFC	114696		20-30 mmHg
XLFC	114629		20-30 mmHg
S	114635		30-40 mmHg
M	114636		30-40 mmHg
L	114637		30-40 mmHg
XL	114638		30-40 mmHg
LFC	114697		30-40 mmHg
XLFC	114639		30-40 mmHg

Fiber contents:

Knee High With Silicone Dot Band

20-30 mmHg: 74% Nylon, 20% Spandex, 6% Silicone
30-40 mmHg: 68% Nylon, 27% Spandex, 5% Silicone

Knee High

15-20 mmHg: 80% Nylon, 20% Spandex

20-30 mmHg: 77% Nylon, 23% Spandex, KH CT; 75% Nylon, 25% Spandex, KH OT

30-40 mmHg: 72% Nylon, 28% Spandex, KH CT; 70% Nylon, 30% Spandex, KH OT

Thigh High ordering information:

Thigh High with Silicone (Closed Toe)



Size	Catalog Number		Compression
	Beige	Black	
S	114822	114826	15-20 mmHg
M	114823	114827	15-20 mmHg
L	114824	114828	15-20 mmHg
XL	114825	114829	15-20 mmHg
S	114208	114212	20-30 mmHg
M	114209	114213	20-30 mmHg
L	114210	114214	20-30 mmHg
XL	114211	114215	20-30 mmHg
S	114216	114220	30-40 mmHg
M	114217	114221	30-40 mmHg
L	114218	114222	30-40 mmHg
XL	114219	114223	30-40 mmHg

Thigh High with Silicone (Open Toe)



Size	Catalog Number		Compression
	Beige	Black	
S	114818		15-20 mmHg
M	114819		15-20 mmHg
L	114820		15-20 mmHg
XL	114821		15-20 mmHg
S	114200		20-30 mmHg
M	114201		20-30 mmHg
L	114202		20-30 mmHg
XL	114203		20-30 mmHg
S	114204		30-40 mmHg
M	114205		30-40 mmHg
L	114206		30-40 mmHg
XL	114207		30-40 mmHg

Thigh High without Silicone (Closed Toe)

Size	Catalog Number		Compression
	Beige	Black	
S	114640		20-30 mmHg
M	114641		20-30 mmHg
L	114642		20-30 mmHg
XL	114643		20-30 mmHg
S	114648		30-40 mmHg
M	114649		30-40 mmHg
L	114650		30-40 mmHg
XL	114651		30-40 mmHg

Thigh High without Silicone (Open Toe)

Size	Catalog Number		Compression
	Beige	Black	
S	114644		20-30 mmHg
M	114645		20-30 mmHg
L	114646		20-30 mmHg
XL	114647		20-30 mmHg
S	114652		30-40 mmHg
M	114653		30-40 mmHg
L	114654		30-40 mmHg
XL	114655		30-40 mmHg

Fiber contents:

Thigh High With Silicone Dot Band

15-20 mmHg: 74% Nylon, 18% Spandex, 8% Silicone

20-30 mmHg: 70% Nylon, 23% Spandex, 7% Silicone

30-40 mmHg: 66% Nylon, 28% Spandex, 6% Silicone

Thigh High Without Silicone Dot Band

20-30 mmHg: 74% Nylon, 26% Spandex, TH CT; 73% Nylon, 27% Spandex, TH OT

30-40 mmHg: 67% Nylon, 33% Spandex, TH CT; 68% Nylon, 32% Spandex, TH OT