



ABSOLUTE SUPPORT *for Him*

Absolute Support for Men includes compression socks that are designed for men who are on their feet a great deal, frequently travel by air, are overweight, have swollen ankles and legs, or have varicose veins.



Men's Microfiber Support Socks


Our most comfortable high compression men's sock at an amazing price. The microfiber provides an exceptional, comfortable fit with added softness in an incredibly durable sock.

10% Lycra Spandex, 90% Soft, Microfiber Nylon

Item # A2011
Medium Support
15-20 mmHg
Closed toe
\$12.95
(6+ \$12.49/ea)

Item # A2012
Firm Support
20-30 mmHg
Closed toe
\$16.95
(6+ \$14.95/ea)

Item # A305
X-Firm Support
30-40 mmHg
Closed toe
\$22.95
(6+ \$19.95/ea)




Men's Rayon Premium Light Support Socks

Our most comfortable 8-15mmHg dress sock ever, with the natural comfort properties of Rayon fiber. With the feel of silk and the texture of cashmere, Ames Walker is proud to introduce our newest, unique addition. The Men's Premium Rayon Support sock, with 8-15mmHg compressions, promises extraordinary comfort and style.

70% Rayon, 25% Nylon, 5% Lycra Spandex

Item # A504
Light Support • 8-15 mmHg • Closed toe
\$11.49 (6+\$10.95/ea)



Cotton Knee Hi Socks Medium Support



Healthier, better-feeling legs begin at the ankles. Graduated support applies pressure and stimulates circulation at the ankle, slowly decreasing as it travels up the leg. Moderate graduated support gently massages and enhances upward circulation in the legs. By far our most popular sock, Absolute Support Knee Highs are clinically proven and physician recommended. Natural cotton construction provides a refreshing natural comfort. Thicker than our regular nylon/lycra dress socks, these socks are warmer in winter, cooler in summer, and ideal as an athletic sock for sports or walking. Reinforced Heel and Toe.

80% Cotton, 13% Nylon, 7% Lycra Spandex



Item # A105
Medium Support • 15-20 mmHg • Closed toe
\$14.95 (6+ \$13.95/ea)

Item # A2105
Firm Support • 20-30 mmHg • Closed toe
\$14.95 (6+ \$13.95/ea)





Microfiber/Cotton Dress Socks Medium Support

Super comfortable dress support socks with graduated support Using modern knitting technology, we have added a layer of cotton inside, next to the skin, providing a softer cooler feel. Soft microfiber nylon stretches to fit comfortably. Moisture-wicking fibers keep feet dry and healthy. We've also added a design to the sock to add a bit of elegance to this very unique product.

80% Microfiber Nylon, 5% Dupont Lycra Spandex, 15% Cotton

Item # A1013
 Medium Support • 15-20 mmHg • Closed toe
\$12.95 (6+\$11.95/ea)

Black	Navy	Khaki	M	L	XL
			7 1/2-10	10 1/2-12 1/2	13 & over

Men's Dress Socks Firm Support

Over-the-calf style provides firm graduated compression. Thinner dress sock with reinforced heel and toe, reciprocated heel pocket and a large toe box for a comfortable fit. Non-constricting comfort top offers great stay up power.

90% Nylon, 10% Lycra Spandex

Item # A104
 Extra Firm Support • 20-30 mmHg • Sandle foot design
\$12.95 (6+\$11.95/ea)

Black	Brown	Navy	White	Khaki	Grey	S	M	L	XL
						7 & under	7 1/2-10	10 1/2-12	12 & over



Men's Firm Support Thigh High With Grip Top

Comfortable, gradient compression designed specifically for a man!

These new thigh highs look and feel like the leading name brands but are significantly discounted! Ribbed design, longer length, flat toe seams and a roomy toe section. Reinforced toe and heel. Easy to wear with a comfortable silicone top band that keeps these thigh highs up without binding.

Item # A2017
 Firm Support • 15-20 mmHg • Closed toe
\$16.95 (4+\$15.95/ea)

M	L	XL
7 1/2-10	10 1/2-12 1/2	13 & over

Black



Men's Firm Support Leotard

Masculine, extra firm support, full length garment for treatment of upper and lower leg circulatory problems. This unique product features a brief type top, with a built in fly opening for men. FINALLY!

Features a two-way stretch design for extra comfort and to make them easier to get on.

80% Nylon, 20% Lycra

Item # A234
 Firm Support • 20-30 mmHg • Closed toe
\$39.95

Black		S	M	L	XL	2XL
	Ankle	8-9 1/2	8 3/4-10 1/2	9 1/2-11 1/2	10-12	11-13
	Calf	12 1/2-14 3/4	14-16	15-17 1/2	16-18 3/4	19 1/2-22 1/2
	Thigh	17-19 1/2	18 1/2-21	19 1/2-22 1/2	21-24 1/2	25 1/2-30
	Height	to 5'8"	to 5'11"	to 6'2"	to 6'4"	to 6'5"

